



Ace Medical Supply, LLC

Right Equipment! Better Life!

Specializing in Manual Wheelchairs, Hospital Beds, and Low Air Loss(LAL) Mattress



Semi- Automatic **Hospital Bed**

(with) Pressure Relieving/Reducing Mattress, &
Full size Bed Rails

At no additional cost for all MEDICARE patients

The Medline Hospital Bed is Lighter, Easier, Stronger and Safer!

Medline Hospital Bed is Lighter

The Medline Full Electric Hospital Bed weighs only 141 lbs. Smooth rolling, 3" casters make the Medline Bed very easy to move. Two of the Medline Hospital Bed casters have top activated brakes to lock the bed in place for added security and peace of mind.

Medline Hospital Bed is Easier

Medline Hospital Bed electronics and motors are washable with standard water pressure. Plastic Medline Bed end panels are easy to wipe down and keep clean. Medline Beds can be easily transported without removal of the motors.

New Medline Bed is Stronger

Medline Beds have a 450 lb. weight capacity with an extra steel plate welded to the middle of the frame. Plastic Medline Bed end panels resist breakage while the heavy duty Medline Hospital Bed remote resists breakage and moisture. Medline Bed frame is painted with three coats of enamel and clear coated to provide a strong finish that is chip and scratch resistant.

Medline Beds are Safer

Medline Beds feature quiet DC motors that keep patients and caregivers safe with a low voltage system that is also washable. Medline Hospital Bed headspring supports are moved out of the way to resist pinching and injuries. Medline Bed foot assembly hinges are also placed inside the frame to significantly reduce any chance of pinching. The Medline Full Electric Hospital Bed includes a lifetime warranty on frame and welds as well as any product defect for the original consumer.

- Bed Width: 34 1/2"
- Bed Height: 14 1/2" min. - 22 1/2" max. (Measurement from floor to bed deck.)
- Dimensions: Overall: 35 1/4" W x 89" L
- Product Weight: 141 lbs.
- Product Weight Capacity: 450 lbs.
- Manufacturer's Warranty: 1 year

***** Medicare does not pay for the Pressure Relieving/Reducing Foam Mattress. However, Ace Medical Supply LLC provides all patients with fully automatic hospital beds even if Medicare do not cover the items.**



MedLite Beds



Lighter,
Easier,
Stronger,
Safer.



MedLite Beds and Accessories

Our Beds are Lighter

- Total semi-electric bed weighs only 137 lbs.
- Smooth-rolling, 3" casters make bed easy to move. Two casters have top-activated brakes to lock in place.
- Optional bed cart is available, and makes storage and transport of the bed a breeze. (MDSBEDCART)

Easier

- Electronics can be washed down with standard water pressure. (No high pressure washers).
- Plastic end panels are easy to wipe down and keep clean.
- Bed can be easily transported without removal of motors.

Stronger

- 450 lb. weight capacity.
- Extra steel plate welded to middle of frame.
- Plastic end panels resist breakage.
- Heavy-duty remote resists breakage and moisture.
- Extended five-year warranty on all electrical and mechanical parts. (Lifetime warranty on frame).
- Frame is painted with three coats of enamel and clear-coated to provide a strong finish that is chip and scratch resistant.

Safer

- Quiet DC Motors keep patients and caregivers safe with a low voltage system.
- Headspring supports moved out of way to resist pinching and injuries.
- Foot assembly hinges placed inside frame to significantly reduce chance of pinching.



- MedLite Beds have:
- Washable Electronics
 - Sealed Durable Panels
 - Low Bed Option


Our Motors are Quieter, Lighter and Washable



Our NEW Remote is Durable and Comfortable



Quality, Extra Steel Plates and Welding Prevent Sagging



Specifications

Dimensions: Length x Width	Spring Deck Width	Spring Deck Height	Section Adjust Head Knee
89" x 35-1/4"	34-1/2"	14-1/2" to 22-1/2"	75° 30°

Bed Frames

Description	Item #	Wt. (lbs.)
Semi-Electric, Two Motor, One-Crank Bed	MDR107002L NEW!	137
Full Electric, Three Motor Bed	MDR107003L NEW!	141
Heavy-Duty Manual, Three Crank Bed (500# Capacity)	MDR107000P	145

